

RADA HYSTEROMETER


The Hysterometer is a sterile, single-use medical device designed to measure the depth of the uterine cavity during gynecological and fertility procedures.



Accurate Uterine Measurement
Precisely Measures Uterine Cavity Depth During Clinical Procedures.



Calibrated Scale
Clear Depth Markings Support Reliable Uterine Assessment.



Gentle Insertion
Smooth Cervical Passage With Minimal Patient Discomfort.



Clinical Guidance
Supports Evaluation And Treatment Planning In Gynecology And Fertility Care.

It allows accurate assessment of uterine length to support clinical evaluation and assist in treatment planning, while enabling smooth cervical passage with minimal patient discomfort.

KEY FEATURES

- Single-use hysteroscope
- Flexible tube with sufficient rigidity for controlled insertion
- Rounded, atraumatic distal tip to minimize cervical and uterine trauma
- Clear length markings up to 16 cm for accurate uterine depth measurement
- Suitable for routine gynecological and fertility procedures
- EO Sterilized

AVAILABLE MODELS

IMAGE	MODEL	DESCRIPTION	TRANSFER TUBE LENGTH	ORDER NO.
	Hysterometer	Standard Straight Model With Graduated Scale	16 cm, graduated	BR344184