

RADA IUI CATHETER



The Intrauterine Insemination (IUI) Catheter is a **sterile, flexible, and safe device** designed for **atraumatic placement of prepared sperm directly into the uterine cavity** during IUI procedures.



Precise Sperm Delivery
Controlled Placement Of Prepared Sperm Into The Uterine Cavity.



Smooth Navigation
Flexible Design Enables Gentle Cervical Passage.



Atraumatic Tip
Soft Rounded Tip Helps Minimize Patient Discomfort.




Reliable Transfer
Optimized Lumen Supports Accurate And Efficient Sperm Transfer.




It allows **smooth navigation through the cervical canal** while minimizing patient discomfort.

Featuring a **soft, rounded distal tip and controlled flexibility**, the catheter ensures **accurate placement under gentle manual guidance, with or without ultrasound assistance**.

Its **optimized lumen and material selection** support efficient sperm transfer while maintaining **structural stability during insertion**.



KEY FEATURES

 Available in three models: Standard, Mandrel, and Curved (Curved Model), with selection based on clinical use and physician preference	 Soft, flexible tube with controlled flexibility for smooth and stable insertion	 Rounded, atraumatic distal tip to reduce uterine and cervical irritation
 Optimized inner lumen for efficient sperm transfer	 Compatible with standard syringes; optional depth markings for guidance	 Suitable for a wide range of patient anatomies
 Adjustable marker for depth control	 Has two side spray holes at the tip	 Ability to connect a lower lock to the hub
 Sterilized with ethylene oxide (EO) and single-use, suitable for fertility clinics and ART laboratories	 Shelf Life 3 Years: Sterility and functionality guaranteed up to three years from manufacture	

AVAILABLE MODELS

IMAGE	MODEL	DESCRIPTION	ORDER NO.
	without mandrel	Standard IUI catheter with soft, flexible tube and rounded distal tip for atraumatic sperm placement. Suitable for routine intrauterine insemination.	BR110134
	with mandrel	IUI catheter with mandrel for added stiffness and control, ensuring precise placement under manual or ultrasound guidance.	BR110133
	Curved Model	IUI catheter with a pre-formed curved tip for enhanced navigation in challenging cervical anatomy, maintaining atraumatic insertion.	BR110135
	With Syringe	IUI catheter supplied with a 1 cc syringe for controlled loading and intrauterine insemination.	